

912-0012 NC-250 IQ/OQ Test Kit

Contents	<p>P/N: 942-0001 1 box of NC-Slide A2 containing 25 glass slides.</p> <p>P/N: 071-0015 1 bottle with 1 mL of fluorescent beads. The solution contains multicoloured 10 µm beads, Tween-20 (0.02%) and NaN3 (2 mM ~ 0.01%) in H2O.</p> <p>CAS no. for Tween-20 is 9005-64-5 EINECS no. for NaN3 is 247-852-1.</p>
Application	<p>NC-250 IQ/OQ Test Kit is part of the <i>NucleoCounter® NC-250™</i> system, and is used for performing installation qualification and operation qualification (IQ/OQ).</p>
Principle	<p>NC-250 IQ/OQ Test Kit contains multicoloured beads and slides that allow testing of the performance of all the channels of the <i>NucleoCounter® NC-250™</i>.</p>
Use	<p>With respect to the detailed procedure, please refer to the appropriate technical notes and user manuals for the <i>NucleoCounter® NC-250™</i> instrument and the <i>NucleoView NC-250 Software™</i>.</p>
Storage	<p>The fluorescent bead solution should be stored at 2-8°C. The <i>NC-Slide A2</i> should be stored at room temperature < 40°C.</p>
Stability	<p>The stability of the kit is 15 months if stored under recommended conditions. Expiry dates are shown on the slide box as well as on the fluorescent bead bottle. The NC-250 IQ/OQ Test Kit is produced 15 months before the expiry date.</p> <p>Discard the NC-250 IQ/OQ Test Kit after performing the OQ procedure.</p>
Safety information	<p>Fluorescent bead solution: In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.</p> <p>If a spill is observed perform a clean-up of the area which may have been in contact with the solution.</p> <p>Use gloves and suitable protective clothing.</p> <p>Please refer to the SDS regarding further safety information.</p>
Disposal	<p>After use, dispose of the NC-250 IQ/OQ Test Kit following national or local laws and regulations, according to the nature of its contents.</p>