
NucleoCounter[®] NC-200[™]

Instrument User Guide

Product No. 991-0203

Revision 1.10

NucleoCounter[®] NC-200[™]

Flexible Cell Counter

Product No. 991-0203 (English)

Revision 1.10

27jan2023

ChemoMetec A/S

Gydevang 43, DK-3450 Allerød, Denmark

Telephone: (+45) 48 13 10 20, Fax: (+45) 48 13 10 21

Internet: www.chemometec.com

E-mail: sales@chemometec.com

Caution!

This equipment must be operated as described in this User Guide and documents referred to herein. Please read the entire guide and referred documents before attempting to use this unit. Please pay attention to that gloves or protective clothing are not worn on the illustrations/pictures shown in this User Guide. However, ChemoMetec A/S does recommend that the user wear suitable protective clothing etc.

Contacting support

Technical information including product literature, answers to questions regarding the operation of the NucleoCounter® NC-200™ not covered in this document and referred documents is available through the following:

- For e-mail support, send questions to NucleoCounter® NC-200™ Technical Support on the address Support@chemometec.com
- Check out the FAQ section under support at www.chemometec.com
- To speak with a Technical Support Specialist, call (+45) 48 13 10 20.

Please note the NucleoCounter® NC-200™ serial number and have it available when contacting ChemoMetec A/S for support. The NucleoCounter® NC-200™ serial number is found on the label affixed to the bottom of the instrument. The version number of the NucleoView NC-200™ software shall also be noted, this can be found on the Help – About menu item in the NucleoView NC-200™ software.

Sales and ordering information

For sales assistance with NucleoCounter® NC-200™ or the NucleoView NC-200™ software, to place an order for a NucleoCounter® NC-200™ or consumables, call (+45) 48 13 10 20, fax (+45) 48 13 10 21, or send e-mail to Sales@chemometec.com

Disclaimer Notices

The material in this document and referred documents is for information only and is subject to change without notice. While reasonable efforts have been made in preparation of this document to provide accurate and timely information, ChemoMetec A/S has no liability and assumes no responsibility or liability resulting from errors or omissions in the content of this document, or from the use of the information contained herein. The information contained in this document is provided on an “as is” basis with no guarantees of completeness, accuracy, adequacy, usefulness, timeliness, or fitness for any specific purpose of the information included or referred to in this document. There may or may not be technical inaccuracies or typographical errors or hazard symbol errors contained in this document, and ChemoMetec A/S will correct such items at its sole discretion as they are discovered.

ChemoMetec A/S reserves the right to make changes in the product design without reservation and without notification to its users.

Copyright Notices

Copyright © ChemoMetec A/S 2011. All rights reserved. No part of this publication and referred documents may be reproduced, stored in a retrieval system, or transmitted in any form or by any means, electronic, mechanical, photocopying, recording or otherwise, without the prior written consent of ChemoMetec A/S, Gydevang 43, DK-3450 Allerød, Denmark.

ChemoMetec and NucleoCounter are registered trademarks owned by ChemoMetec A/S. NucleoCounter, NucleoCassette, Via1-Cassette, Via2-Cassette and NucleoView are trademarks of ChemoMetec A/S.

All other trademarks are the property of their respective owners.

Declaration of Conformity (CE) - *EXAMPLE*

Name of product: NucleoCounter® NC-200™

Type: Flexible Cell Counter

Other identifying data: Product number 900-0201

We declare under our sole responsibility that the products, to which this declaration relates, are in conformity with the Council Directives on the approximation of the laws of the EEC Member States relating to the following:

Directives Electromagnetic Compatibility (EMC) Directive 2014/30/EU
 Restriction of the use of certain hazardous substances (RoHS) Directive 2011/65/EU
 Low Voltage Directive (LVD) 2014/35/EU
 Machinery Directive (MD) 2006/42/EC

Harmonized Standards

EMC: EN 61326-1:2013 Electrical equipment for measurement, control, and laboratory use – EMC requirements – part 1: General requirements
 EN 61000-3-2:2014 Limits for harmonic current emissions
 EN 61000-3-3:2013 Limitation of voltage changes, voltage fluctuations and flicker

RoHS: EN 63000:2018 Technical Documentation for the assessment of electrical and electronic products with respect to the restriction of Hazardous substances

LVD: EN 61010-1:2010+A1:2019 Safety requirements for electrical equipment for measurement, control, and laboratory use - Part 1: General requirements
 EN 61010-1-081 Particular requirements for automatic or semi-automatic laboratory equipment for analysis and other purposes.

MD: EN 12100:2010 Safety of machinery; risk assessment and risk reduction
 EN 13849-1:2015 Safety-related parts of control systems – Part 1: General principles for design

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Date & Signature will appear on the instrument's actual Declaration of Conformity statements, accompanying each new instrument order.

Name and address of manufacturer:

Date: 2022-08-25

ChemoMetec A/S
 Gydevang 43
 DK-3450 Allerød
 DENMARK

Declaration of Conformity (UKCA) - *EXAMPLE*

Name of product: NucleoCounter® NC-200™
Type: Flexible Cell Counter
Other identifying data: Product number 900-0201

We declare under our sole responsibility that the products, to which this declaration relates, are in conformity with the relevant legislation of United Kingdom:

Regulations Electromagnetic Compatibility Regulations 2016
 The Restriction of Use of Certain Hazardous Substances in Electrical and Electronic Equipment Regulations 2012
 The Electrical Equipment (Safety) Regulations 2016
 The Supply of Machinery (Safety) Regulations 2008

Designated Standards

EMC: BS EN 61326-1:2013 Electrical equipment for measurement, control, and Laboratory use – EMC requirements - General requirements
 BS EN 61000-3-2:2014 Limits for harmonic current emissions
 BS EN 61000-3-3:2013 Limitation of voltage changes, voltage fluctuations and flicker

RoHS: BS EN IEC 63000:2018 Technical documentation for the assessment of electrical and electronic products with respect to the restriction of hazardous substances

LV: BS EN 61010-1:2010+A1:2019 Safety requirements for electrical equipment for measurement, control, and laboratory use – General requirements
 BS EN 61010-1-081 Particular requirements for automatic or semiautomatic laboratory equipment for analysis and other purposes

MR: BS EN 12100:2010 Safety of machinery; risk assessment and risk reduction
 BS EN ISO 13849-1:2015 Safety-related parts of control systems – General principles for design

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Date & Signature will appear on the instrument's actual Declaration of Conformity statements, accompanying each new instrument order.

Name and address of manufacturer:

Date: 2022-08-25

ChemoMetec A/S
 Gydevang 43
 DK-3450 Allerød
 DENMARK

WEEE directive information – Europe only


Correct Disposal of This Product

(Waste Electrical & Electronic Equipment) - Europe only

This marking shown on the product or its literature, indicates that it should not be disposed together with other household wastes at the end of its working life. To prevent possible harm to the environment or human health from uncontrolled waste disposal, please separate this from other types of wastes and recycle it responsibly to promote the sustainable reuse of material resources.

Business users should contact their supplier and check the terms and conditions of the purchase contract. This product should not be mixed with other commercial wastes for disposal.

This information is listed in “[Appendix B: WEEE directive information in more EU languages](#)”.

Introduction and intended use

The NucleoCounter® NC-200™ is a compact cell counter intended to be used for precise cell counting within the life science and pharmaceutical industries. The NucleoCounter® NC-200™ system is intended for research use only, and not for diagnostic use. The NucleoCounter® NC-200™ uses disposable cassettes from ChemoMetec.

The NucleoCounter® NC-200™ is a highly advanced cell counter based on fluorescence microscopy using CMOS camera technology with low magnification and advanced image analysis. Different properties of the individual cells can be analyzed with the three integrated light sources (365nm, 500nm and dark field) and with two emission bands. Spatial position is maintained from filter to filter enabling direct comparison of the different parameters detected or quantified on each individual cell. The low magnification (×1.3) guarantees a high analysis volume and thus ensures a high statistical performance.

Specially designed sample illumination ensures homogeneous illumination of the complete measuring area and thereby accurate quantification of fluorescence signal is achieved. Since no tubes, valves, hoses etc. are present the instrument requires little maintenance.

The NucleoCounter® NC-200™ is improved in relation to NucleoCounter® NC-100 and shares some of the strong features with the NucleoCounter® NC-3000. In normal viability analysis the sample volume is as high as 1.4 µl and the instrument is also able to count cells in aggregates.

The NucleoCounter® NC-200™ part number at ChemoMetec is 900-0201.

Warnings and precautions

Whenever the ⚠ symbol appears on the NucleoCounter® NC-200™ instrument, it indicates that the manual must be consulted for precautions and warnings.

Power and cables

Use the shielded USB cable supplied with the NucleoCounter® NC-200™ to ensure that appropriate EMI classification is maintained for the intended environment.

⚠ The USB interface connector of the NucleoCounter® NC-200™ must only be connected to SELV circuits. External computing devices connected to the USB interface connector of the NucleoCounter® NC-200™ must comply with the standards, UL 1950 and IEC/EN 60950.

Table 1 Description of the NucleoCounter® USB interface connector

Pin no.	Name	Maximum Voltage level ¹
1	+5V	+5 VDC
2	D-	+3.5 VDC
3	D+	+3.5 VDC
4	DGND	0 VDC
Metal enclosure	Shield (connected to DGND)	0 VDC

NucleoCounter® NC-200™ is powered by an external 24VDC power supply. For safe use, please follow the instructions for connecting the power supply.

⚠ The NucleoCounter® NC-200™ shall only be used with an external power supply purchased at ChemoMetec.

⚠ The detachable DC power supply cord and appliance inlet of the external power supply are considered the disconnecting device.

⚠ The mains supply cord and plug of the external power supply shall comply with any national regulations.

⚠ The user should be made aware of that, if the NucleoCounter® NC-200™ and the external power supply are used in a manner not specified by the manufacturer, the

¹ In normal operation mode (refers to Pin no. 4)

protection provided by the NucleoCounter® NC-200™ and the external power supply may be impaired.

Electromagnetic interference

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment to an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Caution! Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Cassettes, Reagents, Solutions and Dispensers

With respect to use and handling of cassettes, reagents, solutions and dispensers, please refer to appropriate package inserts for these items.

Caution! When using a bottle-top dispenser: To protect against accidental splashes, protective clothing, eye protection and gloves must be worn when using potentially hazardous liquids.

General

Any biological material should be handled as if it can transmit infectious disease and disposed of with proper precautions according with federal, state and local regulations.

Avoid specimen contact with skin or mucous membranes.

Never pipette by mouth.

Avoid cross contamination of the samples when preparing the samples. This can compromise the quality of the results.

Table of contents

DECLARATION OF CONFORMITY (CE) - <i>EXAMPLE</i>	V
DECLARATION OF CONFORMITY (UKCA) - <i>EXAMPLE</i>	VI
WEEE DIRECTIVE INFORMATION – EUROPE ONLY	VII
INTRODUCTION AND INTENDED USE	VIII
WARNINGS AND PRECAUTIONS	IX
1 INSPECTION, UNPACKING AND INSTALLATION OF EQUIPMENT	3
1.1 OPERATION OF THE NUCLEOCOUNTER® NC-200™	3
2 CELL COUNTING WITH THE NUCLEOCOUNTER® NC-200™	4
2.1 THE NUCLEOCOUNTER® NC-200™ INSTRUMENT	4
2.1.1 <i>Fluorescence microscope</i>	5
2.1.2 <i>The Cassette fixture</i>	5
2.1.3 <i>The RUN Button and LED indicator</i>	5
2.2 THE VIA1-CASSETTE™	6
2.2.1 <i>Loading the Via1-Cassette™</i>	6
2.2.2 <i>Sample volume</i>	7
2.2.3 <i>Cassette handling</i>	7
3 KEY, LED INDICATOR AND PC CONTROL	8
3.1 INTERACTIVE CONTROLS	8
4 START-UP	9
4.1 POWER ON, POWER OFF	9
4.1.1 <i>Starting up</i>	10
5 OPERATION OF THE NUCLEOCOUNTER®	11
5.1 READY MODE	11
5.2 ANALYSIS WITH THE CASSETTES	11
5.2.1 <i>Actuator movements</i>	12
5.3 RESULT PRESENTATION	12
5.4 PROTECTION OF THE OPTICAL SYSTEM	12
6 NUCLEOVIEW NC-200™	13
7 MAINTENANCE OF NUCLEOCOUNTER® NC-200™	15
7.1 CLEANING	15
7.1.1 <i>Instrument casing</i>	15
7.1.2 <i>Cassette fixture</i>	15
7.1.3 <i>Spill of liquid</i>	16
7.2 DECONTAMINATION	16
8 PC- MAINTENANCE, ERROR MESSAGES AND POWER-ON FAILURE	17

8.1	REQUIRED HARD DISK SPACE FOR DATA STORAGE	17
8.2	ERROR MESSAGES.....	17
8.3	POWER-ON FAILURE.....	17
9	TECHNICAL SPECIFICATIONS	19
9.1	THE NUCLEOCOUNTER® NC-200™.....	19
9.2	SAMPLE MEDIAS	19
9.3	STANDARDS	20
10	EQUIPMENT AND ACCESSORIES.....	21
11	GENERAL DISCLAIMER	22
	APPENDIX A: DESCRIPTION OF CHANGES FROM LATEST REVISION.....	23
A.1	NEW:.....	23
A.2	CHANGES:.....	23
A.3	REMOVED:	23
	APPENDIX B: WEEE DIRECTIVE INFORMATION IN MORE EU LANGUAGES	24

1 Inspection, Unpacking and Installation of Equipment

Upon receiving the order from ChemoMetec, the box or boxes should be carefully inspected for any damage that may have occurred during shipping. Any damage must be reported to the carrier and to ChemoMetec immediately.

Unpack the boxes, saving the packing materials for possible later use.

Verify that the ChemoMetec packing list received refers to the correct and ordered materials, and that nothing is missing.

If any part of the order was damaged during shipping, is missing or fails to operate, please contact ChemoMetec immediately.

1.1 Operation of the NucleoCounter® NC-200™.

The NucleoCounter® NC-200™ Installation Qualification (IQ) Application Note (# 994-0205) describes the IQ protocol. The NucleoCounter® NC-200™ Operational Qualification (OQ) / Zero Count Application Note (# 994-0206) describes the OQ/Zero Count protocol. Execution of the IQ protocol will ensure that the instrument is installed correctly. The OQ protocol verifies that the instrument is operating correctly.

Once IQ, OQ and Zero Count protocols are successfully completed, the NucleoCounter® NC-200™ is ready for use.

2 Cell counting with the NucleoCounter® NC-200™

The NucleoCounter® NC-200™ is a bi-color fluorescent cell-counting instrument. The NucleoCounter® has the advantage that only short training time is required before the user can run highly advanced cell counting protocols. The instrument itself is calibration free. The individual calibration parameters for each cassette will be read and used by the instrument during analysis.

The complete list of protocols is shown in the document; Protocol list for the NucleoCounter® NC-200™. This covers protocols such as:

- Count of Aggregated Cells Assay
- Viability and Cell Count – Aggregated Cells Assay
- Viability and Cell Count Assay

New protocols will be developed on an on-going basis and the updated list can be seen at the web site: www.chemometec.com

2.1 The NucleoCounter® NC-200™ instrument

The NucleoCounter® NC-200™ is shown in Figure 1 below.


Figure 1. NucleoCounter® NC-200™

2.1.1 Fluorescence microscope

The NucleoCounter® NC-200™ is a fluorescent microscope with a magnification of approximately x1.3. It is equipped with three light sources (two for fluorescent use and one for dark-field analysis). It has a dual-band emission filter.

At present Via1-Cassette™ or Via2-Cassette™ are used for the NucleoCounter® NC-200™.

Several different protocols can be run on the NucleoCounter® NC-200™. All of them are documented in application notes describing the procedure to be used for the protocol. The basis of all these protocols is to provide the operator with a simple set-up to perform an extremely advanced analysis in a very easy and flexible manner so that all personnel can be trained in a very short time.

The NucleoCounter® NC-200™ runs together with the NucleoView NC-200™ software. See the Help Menu in the NucleoView NC-200™ software and the “Install_Guide” file on the supplied USB Flash drive for the list of PC Specifications.

2.1.2 The Cassette fixture

The cassette fixture holds the cassette during analysis in the instrument. Refer to section “7.1 Cleaning” for guidelines on how to clean this fixture if required. A stepper motor moves the cassette piston after the RUN button has been activated.

2.1.3 The RUN Button and LED indicator

On the front of the NucleoCounter® NC-200™ there is a RUN button. Activating this will run the last performed type of analysis on the instrument. In this way it is easy to perform a number of identical analyses in a routine manner.

On the rear of the instrument, a USB data connection (rev 2.0 or higher), a DC Power Connection and a DC fuse holder can be found. Please take care of these areas and keep them free of liquids and foreign objects at all times.

2.2 The Via1-Cassette™


Figure 2. The Via1-Cassette™. Here fluorescent dyes AO (Acridine Orange) and DAPI (4',6-diamidino-2-phenylindole) are immobilized in the first part of the flow system.

2.2.1 Loading the Via1-Cassette™

The Via1-Cassette™ (Figure 2) is loaded by gently pressing the white piston down to the top of the cassette, which creates a partial vacuum in the flow system. The tip of the cassette must be immersed into a sample mixture when pressing the piston, resulting in the sample being loaded into the flow system. Approximately 60 µl is loaded into the flow system filling the first part of the channel in the cassette. The liquid is transferred to the measurement chamber when the cassette is inserted into the NucleoCounter® NC-200™ instrument and an analysis is initiated. Please note that approximately 100 µl in a normal micro centrifuge tube is recommended to reduce the risk of air being loaded into the system.

Figure 3. Loading of the Cassette

The cassettes to be used with NucleoCounter® NC-200™ are highly specialized devices, which are fundamental to several of the unique features of the NucleoCounter® NC-200™. For instance, in the Via1-Cassette the fluorescent dyes AO (Acridine Orange) and DAPI (4',6-diamidino-2-phenylindole) are immobilized inside the flow channels of the disposable cassette. When the cassette has been loaded with approximately 60 µl of sample mixture, the dyes are dissolved and mixed with the sample, staining the cells.

After being placed in the NucleoCounter® NC-200™, the stained mixture is automatically transferred to the measurement chamber, where the fluorescent image is recorded. In a cassette the measurement volume is as large as 1.4 µl.

After analysis, the sample and dyes are contained inside the cassette, which can be safely discarded. This offers unique and safe sample handling and disposal.

As the Via1-Cassette (and all other cassettes) contains the entire flow system as well as measurement chamber, only basic cleaning and maintenance of the NucleoCounter® NC-200™ instrument itself is needed.

2.2.2 Sample volume

During production, a calibration parameter for every cassette is found, which specifies the precise depth or thickness of the measurement chamber of the device. The calibration parameter is used by the instrument during analysis and the volume analyzed is determined by multiplying the imaged area with the depth of the measurement area. The imaged area is only dependent on the optics of the instrument. Therefore, the area is constant and specific for each instrument.

The calibration parameter is printed on every single cassette as a dot code describing actual thickness of the measurement chamber of the cassette.

2.2.3 Cassette handling

The image analysis method of the NucleoView NC-200™ offers considerable stability, when eliminating non-cellular objects from the image, including scratches and smears on the cassette windows, and thus producing valid results even under extreme conditions. Nevertheless, the quality of the cell count is best ensured by keeping the windows of the cassette as clean as possible.

Caution! Be careful when attempting to wipe off the surface of the cassette to remove any foreign object, since the plastic material of the cassette can be scratched. In section 5.2 “Analysis with the cassettes” see how to load and remove the cassette.

3 Key, LED indicator and PC control

3.1 Interactive Controls


Figure 4. The front of the NucleoCounter® NC-200™

Figure 4 above shows the front of the NucleoCounter® NC-200™. Two operational areas are present:

- Lid with cassette fixture behind
- Run button and a Tri color status LED indicator

The lid is located in the first area where the sample, contained in a cassette is placed. When the lid is opened, a cassette can be placed in the fixture.

The other area is where the Run button and LED indicator are placed. Once a protocol has been selected in the software (see detailed description in the NucleoView NC-200™ Help), an analysis using this protocol can be performed by activating the Run button.

The LED indicator shows the status of the instrument.

- Green - the NucleoCounter® NC-200™ is ready for next operation.
- Red - the NucleoCounter® NC-200™ is busy.
- Orange - the NucleoCounter® NC-200™ is not in connection with the NucleoView NC-200™ software.

4 Start-up

4.1 Power on, power off

On the rear of the instrument, two sockets and one fuse holder are located. When seen from the rear the USB 2.0 port is to the left and the fuse holder and DC Power socket is at the right. The NucleoCounter® NC-200™ is powered by a 24VDC external power supply. To connect the power supply to the NucleoCounter®, connect it to the DC Power Plug (see Figure 5), and plug the mains power cord of the power supply into a wall outlet.

Figure 5. Rear side of the NucleoCounter® NC-200™

The NucleoCounter® NC-200™ is turned "On" once the wall outlet is "On". To power it "Off", either remove the DC Power Plug or turn the wall outlet "Off".

Figure 6. USB 2.0

Figure 7. DC Power Socket

4.1.1 Starting up

The NucleoCounter® NC-200™ is connected to a PC by inserting the respective plugs of a USB cable into the instrument and the PC.

Once the NucleoView NC-200™ program is launched the stepping motor for the piston actuator performs a homing procedure to ensure correct positioning of the actuator.

After this initialization, the status indicator on the front will turn green and the NucleoCounter® NC-200™ is ready for operation.

5 Operation of the NucleoCounter®

5.1 Ready mode

Once the NucleoCounter® NC-200™ and the NucleoView NC-200™ software are ready, analysis can begin. In the NucleoView NC-200™ software (Figure 8), the browse menu can be used to select the protocol to be performed. Once selected, the run button on the front of the instrument (see Figure 4) or the run button in the software can be used to execute the analysis.


Figure 8. NucleoView NC-200™

5.2 Analysis with the cassettes

To analyze with a cassette:

1. Load the cassette with an appropriate sample according to the protocol selected.
2. Place the cassette into the cassette fixture and close the lid (Figure 9).
3. Press the Run² button.

After analysis, the cassette is removed from the fixture and disposed of as biological waste according to national or regional laws or regulations regarding its contents.

² Before this action, the correct protocol has to be selected in the NucleoView NC-200™ software.

Figure 9. Illustration of how to hold and insert the cassette in the fixture

5.2.1 Actuator movements

When a cassette is analyzed after Run is pressed, the piston actuator will move the piston rod down through the cylinder of the cassette. In this way, the sample material is transported to the cassette measurement chamber. The form and dimensions of the flow channels facilitate effective mixing of the sample. Sensors monitor sample flow. When the sample has been correctly loaded into the measurement chamber, the piston is stopped.

5.3 Result presentation

Depending on the type of analysis performed the result will be presented on the PC within a few seconds and the fluorescent image(s) will be stored on the hard disk. Typical process time is less than a minute for viability analysis using one Via-1 Cassette.

5.4 Protection of the optical system

The lid covers the fixture. This lid is designed to protect the area from dust and other potential contaminants, and to keep external light from interfering with the fluorescent image recording. The lid is mounted with a magnetic hinge.

Unless inserting or removing a cassette or cleaning the fixture, the lid should be closed at all times to protect the optical system of the NucleoCounter® NC-200™.

6 NucleoView NC-200™


Figure 10. NucleoView NC-200™ with a typical result after analysis

The NucleoView NC-200™ software controls the NucleoCounter® NC-200™ instrument and image analysis of the recorded images. Several options for data presentation are available in the NucleoView NC-200™ software (Figures 10 and 11). The main results of the analysis are presented directly on the user interface of the software. An extended result table is available in the result tab. Choosing files in the browse window directly on the user interface will display the results of the selected file.


Figure 11. A typical analysis with overlay function enabled

The events registered as cells can (when enabled in the protocol) be inspected by launching the Protocol Adaption Editor in the Tools menu and selecting to show events included in the count prior to running an analysis. This will launch a plot manager that shows a scatter plot with a gate that encloses the events recorded as cells (Figure 12). Furthermore, two histograms will be displayed showing the intensity and area of the events in the picture.


Figure 12. The Plot Manager showing a scatter plot and two histograms

A USB Flash drive supplied with the instrument contains the software, drivers, instrument specific data etc. On the same USB Flash drive, you can also find the NucleoView NC-200™ User Guide. Run the file "Install_Guide" and follow the directions to the NucleoView NC-200™ User Guide. Once the program is installed the NucleoView NC-200™ User Guide can be displayed by following the directions: Menu – Help – Software Users Guide.

Refer to the NucleoView NC-200™ User Guide for additional information.

7 Maintenance of NucleoCounter® NC-200™

7.1 Cleaning

7.1.1 Instrument casing

Depending on the environment in which the NucleoCounter® NC-200™ is operated, it is suggested that regular cleaning of the casing is carried out. When cleaning the casing, it is recommended to use a soft moist cloth and gently wipe the surface. Any contamination, which does not come off immediately, should be rubbed gently with a cloth wetted with mild detergent. Never use organic solvents or aggressive detergents to clean the exterior of the NucleoCounter® NC-200™ as this might damage the surface. Daily cleaning of the casing can be performed by using a mild detergent.

7.1.2 Cassette fixture

When cleaning the cassette fixture, great care must be taken against introducing any liquid or dust into the NucleoCounter® NC-200™. Any liquid that enters the interior of the NucleoCounter® NC-200™ can damage the optical parts and thus compromise the quality of the cell counts.

The black surface of the cassette fixture is anodized aluminum and can therefore be wiped with a clean, dry and dust-free cloth.

An object on the surface of an optical component can influence the collected image. A contaminant will normally be visible as a faint object in the image. In case the contamination is on the cassette it will not be visible when a new cassette is analyzed. If the contamination is present on every image taken, it is most likely that other optical parts are contaminated. In this case contact ChemoMetec. An example of a relatively large contamination is given in Figure 13.

Figure 13. An image showing a contamination of the optical system, visible as a white cloudy phenomenon in the upper right-hand corner.

The presence of a contaminant (e.g., dust particles) will normally not influence the counting of cells. The NucleoCounter® NC-200™ will distinguish between cells and contaminants, since a cell generally is significantly smaller than contaminants.

7.1.3 Spill of liquid

If liquid has been spilled on the instrument it can contaminate elements of the optical system. This contamination can show up in the image in several ways and have many causes, but a common feature is that it cannot be removed by the use of compressed air.

Even though such phenomena are clearly visible on the image it only rarely affects the results of the NucleoCounter® NC-200™. Since it is possible to damage the optical system while attempting to clean the system, it is recommended that ChemoMetec should be consulted before the user attempts such cleaning.

7.2 Decontamination

In case decontamination of the NucleoCounter® NC-200™ is required, please contact ChemoMetec before taking action.

8 PC- maintenance, Error messages and Power-on failure

8.1 Required hard disk space for data storage

The CM files stored on the hard disk can vary from 6MB to 30MB. Therefore, it is important from time to time to check whether enough free space is present on the hard disk. ChemoMetec recommends a minimum of 10GBs free space on the hard disk.

It is recommended to make a periodic archiving of CM files to external storage. Those archived CM files can then be removed from the hard disk of the PC.

The periodic archiving can be made on either the creation date of the CM file, or when the free disk space is below a certain amount, or when the disk space used for CM files exceeds a certain amount. Contact the local IT Administrator in order to get information about procedures for external data storing.

A good example of the archiving procedure could be based upon:

- CM files older than three months are archived.
- CM files uses more than 50GBs of disk space, archive the oldest CM files until less than 20GB is used.
- When less than 10 GBs free space, archive the oldest CM files until less than 20GB is used.

Refer to the NucleoView NC-200™ User Guide for additional information.

8.2 Error messages

Under certain conditions the NucleoView NC-200™ software will display error messages during operation. Corrective actions can be found in the Help section of the NucleoView NC-200™ software, but if they do not correct the errors, contact ChemoMetec or the local distributor.

8.3 Power-on failure

If the NucleoCounter® NC-200™ cannot be turned on, please verify that the power supply is connected to the NucleoCounter® NC-200™ and a working mains outlet.

If the NucleoCounter® NC-200™ still does not turn on, inspect and/or replace the fuse as described below.

The fuse holder is located at the rear side above the power plug, see Figure 14.

Figure 14. Location and removal of the fuse

Disconnect the power supply from the NucleoCounter®. Then remove the fuse holder using a screwdriver as indicated in Figure 14 (turn counter clockwise). Inspect the fuse. If it needs replacement use one as described by ChemoMetec.

If the NucleoCounter® NC-200™ does not turn on after replacement of the fuse, or if the fuse was not broken, please contact ChemoMetec or the local distributor for further instructions.

9 Technical specifications

9.1 The NucleoCounter® NC-200™

Optics	Lens with $\times 1.3$ magnification, 1/2" CMOS
Sample device	Via1-Cassette™
Excitation (nm)	Two LED light sources with peaks at 365 nm, 505 nm and dark field light
Emission (nm)	A single dual-band emission filter: 410-460 nm and 540-650 nm
Sample consumption	60 μ l
Analyzed volume	1.4 μ l sample
Optimal Range	5×10^4 - 5×10^6 cells/ml
Cell types	Mammalian cells
PC Platform	Microsoft Windows 10 in 32 or 64 bits. At least one USB 2.0 port must be available. Screen minimum 1024 $\times$ 768. At least 2 GB RAM and 10 GB free disc space.
Analysis time	Typical 50-60 seconds
Data presentation	Images, tables, histograms, scatter plot
Weight and Dimensions	4.5 kg, H 23 cm, W 40 cm, D 25 cm
Supply Voltage	24 VDC (100-240V~ 50-60Hz)
Power Consumption	7/20W (Ready mode/peak)
Operation and storage conditions	Maximum relative humidity 80 percent for temperatures up to 31 °C decreasing linearly to 65 percent relative humidity at maximum 35 °C; minimum temperature 15 °C.

9.2 Sample medias

Storage and stability	Refer to package insert for the appropriate sample media type.
-----------------------	--

9.3 Standards

Please refer to the section [Declaration of Conformity \(CE\) - EXAMPLE](#) or [Declaration of Conformity \(UKCA\) - EXAMPLE](#) for an overview of compliance to the EU and UK standards.

10 Equipment and Accessories

On www.chemometec.com an updated list of Equipment and Accessories for the NucleoCounter® NC-200™ is available.

11 General Disclaimer

The information contained herein is to the best of our knowledge accurate and complete. However, cell species and cell environments may vary in property. Therefore, systematic and/or random deviation between estimates obtained by the NucleoCounter® and other cell counting methods may occur. As such, nothing contained or stated herein including results obtained from use of NucleoCounter® NC-202™ or Via2-Cassette™ shall be construed to imply any warranty or guarantee. ChemoMetec and affiliated companies shall not be held liable for damages and customers shall indemnify ChemoMetec and affiliated companies against liability resulting from use of potentially inaccurate data generated by the NucleoCounter®. It is recommended that all results obtained using the NucleoCounter® NC-202™ are validated against appropriate reference methods and/or traditional laboratory methods at regular intervals.

Appendix A: Description of changes from latest revision

This chapter describes changes from latest revision of this User guide.

A.1 New:

- Addition of Declaration of Conformity (UKCA) – *EXAMPLE*.

A.2 Changes:

- Disclaimer Notices updated.
- Declaration of Conformity (CE) – *EXAMPLE*
 - Updated to comply with current directives.
- Technical Specifications, 'Dimensions' has been modified for accuracy.
- Technical Specifications, 'Analysis time' has been modified for accuracy.
- Section 1 and 1.1 updated.
- Section 7.3 updated.

A.3 Removed:

- N/A

Appendix B: WEEE directive information in more EU languages

United Kingdom

Correct Disposal of This Product (Waste Electrical & Electronic Equipment) - Europe only

This marking shown on the product or its literature, indicates that it should not be disposed together with other household wastes at the end of its working life. To prevent possible harm to the environment or human health from uncontrolled waste disposal, please separate this from other types of wastes and recycle it responsibly to promote the sustainable reuse of material resources.

Business users should contact their supplier and check the terms and conditions of the purchase contract. This product should not be mixed with other commercial wastes for disposal.

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