

910-3008 Solution 8 - DAPI•PBS**25 ml**

Contents	<p>Solution 8 - DAPI•PBS (25 ml).</p> <p>The solution contains 4',6-diamidino-2-phenylindole (DAPI; 1 µg/ml ~ 3 µM ~ 0.0001%) in phosphate buffered saline (PBS; pH 7.4) supplemented with sodium azide (NaN₃; 0.1 mg/ml ~ 1.5 mM ~ 0.01 %).</p> <p>The EINECS no. for DAPI is 249-186-7.</p> <p>The EINECS no. for NaN₃ is 247-852-1.</p>
Application	<p>Solution 8 - DAPI•PBS stains the nucleus of necrotic and late-apoptotic cells.</p>
Principle	<p>Solution 8 - DAPI•PBS stains the nucleus of cells with compromised plasma membranes allowing enumeration of non-viable cells.</p>
Use	<p>The solution is for research use only (RUO) and is not for diagnostic or therapeutic uses.</p>
Storage	<p>Store at 2-8°C.</p>
Stability	<p>Unopened bottles can be stored until the expiry date indicated on the bottle. The solution is produced 15 months prior to the expiry date.</p> <p>We recommend using the solution no later than the expiry date or two (2) months after opening the bottle, whichever comes first.</p>
Safety information	<p>In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. If a spill is observed, perform a clean-up of the area which may have been in contact with the solution.</p> <p>Use gloves and suitable protective clothing.</p> <p>Please also refer to the safety data sheet (SDS) for safety information.</p>
Disposal	<p>After use, dispose of the bottle following national or local laws and regulations, according to the nature of its contents.</p>